

# Barrier-Free® Examination Tables

## Comparison Chart

Features	Ritter 222	Ritter 223	Midmark 622	Midmark 623	Midmark 625
Power seat height	•	•	•	•	•
Power back rest	Manual	•	Manual	•	•
Power foot support	Manual	Manual	Manual	Manual	Manual
UL 60601-1 approved	•	•	•	•	•
CAN/CSA C22.2 60601.1-M90 approved	•	•	•	•	•
Duplex Hospital-grade Electrical Outlets	•	•	•	•	•
Hideaway stirrups and treatment pan	•	•	•	•	•
Exam Assistant™ Front Drawer	•	•	•	•	•
Paper Roll Holder	•	•	•	•	•
Paper Roll Storage (up to 4 rolls)	•	•	•	•	•
Rear Storage Pod	Optional	Optional	Optional	•	Optional
Pelvic Tilt/Drawer Heater	Optional	Optional	Optional	Optional	Optional
Steel-reinforced Back/Seat section	•	•	•	•	•
Foot control	•	•	•	•	•
Hand control		Optional		•	•
Auto-Return				•	•
Quick Exam					•
Motion Profiling					•
Seamless, removable upholstery	•	•	•	•	•
Ultra-premium, removable upholstery			Optional	Optional	Optional
Premium Comfort System™			•	•	•
Heated Upholstery			Optional	Optional	Optional
Upholstery Color Options	8	8	27	27	27
<b>Weights &amp; Dimensions</b>					
Patient Weight Capacity	400 lbs.	400 lbs.	400 lbs.	400 lbs.	650 lbs.
Minimum Seat Height	18"	18"	18"	18"	18"
Maximum Seat Height	37"	37"	37"	37"	37"
Minimum Seat Height in FLAT-back position	18"	18"	18"	18"	18"
Manual Pelvic Tilt Angle	7°	7°	7°	7°	7°
Back Upholstery Width	28.5"	28.5"	27"	27"	27"
Seat Upholstery Width	28.5"	28.5"	28.5"	28.5"	28.5"
Upholstery Length	59" - 75"	59" - 75"	59" - 75.5"	59" - 75.5"	59" - 75.5"
Footrest Upholstery Length	16"	16"	16.5"	16.5"	16.5"
Footrest Upholstery Width	19.75"	19.75"	19.75"	19.75"	19.75"

### Optional Accessories Available

Adjustable Arm System, Assist Handles, Supply Assistant, Rear Storage Pod, Articulating & Standard Knee crutches, Procto rest, Urology drain pan, Welch Allyn light bracket, PC work surface, Stainless steel treatment pan, Flat upholstery, Thick footrest.

For product inquiries, please contact Midmark Medical Customer Service at 1-800-MIDMARK, extension 88921.  
 For more information or a demonstration, contact Midmark, 60 Vista Drive, P.O. Box 286, Versailles, Ohio 45380-0286.  
 Call: 1-800-MIDMARK Fax: 1-800-365-8631 midmark.com

